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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

## Hip orthosis

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/ or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields. Wearing a knee-length hose under the Kafo module and a sports shoe with laces is recommended. Do not apply in direct contact with open wounds.

### WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.

To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician.

The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation.

The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic device is only guaranteed when all the components are properly in use.

### MAINTENANCE

- |               |                      |
|---------------|----------------------|
| Do not bleach | No chemical cleaning |
| Do not iron   | Do not tumble-dry    |

### Washing instructions:

- Padding: Wash by hand in lukewarm water (max. 30°C) with neutral soap; rinse thoroughly. Dry away from heat sources
- Rigid parts: Clean with a sponge soaked in lukewarm water (max. 30°C) and neutral soap. Dry with a cloth
- Check at regular intervals to make sure that the screws for flexion/ extension and the abduction setting are tightened
- Do not dispose of the product or any of its components into the environment

### INDICATIONS

- Pre and/or post-surgery (arthroplasty or hip reison)
- Inoperable or insufficiently consolidated femoral neck fractures
- Abduction and flexion-extension control for patients with potential joint problems

### CONTRAINDICATIONS

Currently no known

### CHARACTERISTICS AND MATERIALS

#### Hipo

- Anatomically shaped and perfectly adaptable plastic parts
- Modular system that permits replacement and adaptation of parts
- Available in 4 sizes, right and left Kinesio hinges
- Abduction/adduction adjustable in steps of 7,5°
- Flexion-Extension settings from 0° to 120° in steps of 15°
- The Kinesio Offset articulations are manufactured at a special bending angle to provide more room between the articulation and the hip.
- REF. 2306 only: Using the dual hinge the position of the thigh component can be corrected medially or laterally

### ADAPTING

#### ADAPTING FOR THE DOCTOR / TECHNICIAN

The hip brace is made up of: Hip holder, Kinesio joint brace and thigh holder.

#### 1 Putting on the hip holder:

- Loosen the four screws of the lumbar clips of the hip holder (but do not remove fully).
- Open the front buckle of the hip holder.
- Position the hip holder in such a way that the bracket for the joint brace is placed along the centre line of the patient's hip (fig. A), and then repeat the same process for the hip holder on the opposite side.
- Check to ensure that the hip holder rests on the iliac crest and the brace fits closely to the waist.
- Exert a light pressure on the brace so that it is pushed close to the hip and tighten the four screws of the lumbar clip (this process is best carried out by two people) (fig. B).
- Check to make sure that the connecting strap is exactly centred (this can be carried out using the notches in the connecting strap).
- Set the stomach strap to the appropriate length and hook into place (fig. C). If necessary, cut off any protruding parts and carefully seal the edge using a small flame.

#### 2 Setting the Kinesio joint brace:

- Remove the clear protective cap after undoing the appropriate fixing screw (fig. D).
- Set the flexion-extension as prescribed by the doctor. Undo the screws, set the degree of flexion-extension and tighten the screws again (fig. E).
- Replace the clear protective cap again and secure using the appropriate screw.
- Loosen the screws for the abduction/adduction joint far enough (but do not remove) that both parts of the brace can rotate freely (fig. F).

#### 3 Attaching the joint brace to the hip holder (fig. G):

- Attach using the supplied screws, do not tighten the screws.
- Align the flexion-extension joint in such a way that the centre point is positioned at the head of the femur neck (to align exactly use the Trochanter major as a reference point). Then tighten both fixing screws on the hip holder.

#### 4 Attaching the thigh holder:

- The bracket for the joint brace must be on the side.
- Ensure that the thigh holder is in a medial position on the condyle and the gap proximal to the bar is at least three to four fingers.
- Close the thigh holder using the appropriate Orthèse pour hanche straps. Before doing so, adjust to the correct length (fig. H). If necessary, cut off any protruding parts and carefully seal the edge using a small flame (cigarette lighter).
- Align the bracket for the joint brace in accordance with the joint brace. Secure the brace at the thigh holder using the two supplied screws. Tighten the screws (fig. I).

#### 5 Setting the adduction/abduction (fig. I):

- Bring the extremity to the abduction or adduction position as prescribed by the doctor.
- In order to secure the required position, tighten the screw for the adduction/abduction joint.

**CAUTION:** If the brace is fitted to a bedridden patient, we recommend setting the hip holder as follows:

- Measure the length of the rear waist circumference (between the centre lines of the hips);
- Loosen the four screws of the lumbar clips of the hip holder (but do not remove fully). Set the hip holder in such a way that the gap (measured on the exterior circumference) between the bracket for the joint brace and the centre point of the opposite pelvic half is equal to half of the waist circumference of the patient, increased by 3-4 cm. Tighten the screws of the lumbar clips. Use the notches on the connecting strap to centre it;
- The angling of the lateral cups of the hip holder opposite the connecting strap must be determined by eye until the patient is able to stand up again.

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