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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM**

Brace for acromioclavicular joint subluxation

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin.
If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully.
We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient.
To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. S1063
Size	universal

Colour: blue - specify right or left

MAINTENANCE
 Do not bleach No chemical cleaning
 Do not iron Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Conservative treatment for stabilization of the acromioclavicular joint after traumatic dislocations of first-second-third degrees (at the discretion of the physician)
- Conservative treatment after fracture of far distal clavicle

CONTRAINDICATIONS

Currently no known

FEATURES AND MATERIALS

- System designed to push clavicle downward and the humerus upward to align acromion and clavicle (exclusive)
- Anatomical, flexible clavicle brace with polyethylene foam padding and rotating joints
- Lateral strap with buckle (orange) for a rapid removal
- Adjustable, comfortable elbow and forearm support, polyethylene foam semi-rigid padding for an easy application
- Adjustable fixing belt snap (white) with Velcro® fastening system
- Rear straps with both constant and pre-set tensioning
- Velour straps, adaptable size

PUTTING ON THE APPLIANCE ADAPTING FOR THE DOCTOR/TECHNICIAN

- 1 Ensure that all straps have been adjusted to maximum length; only open the white clip buckle closure (fig. A) and the orange-coloured hook closure (strap number [3]) (fig. B).
- 2 Position the elbow pocket with forearm support at elbow height (make sure that it is not too low). The elbow pocket must be on the outside with the Velcro® closure facing upwards.
- 3 Close the waist strap with the white clip buckle closure (fig. C) and tighten using the appropriate Velcro® closure (fig. D). If the strap is too long, remove the Y-shaped Velcro® closure, cut the strap to length and replace the Y-shaped Velcro® closure again.
- 4 Align the elbow pocket with the luxated clavicle so that they form a line and insert the elbow into the elbow pocket (fig. E). Regulate the size of the elbow pocket using the appropriate Velcro® closure (fig. F).
- 5 The wrist strap (number [2]) is attached around the forearm and closed using the Velcro® closure at the wrist (fig. G). The wrist strap can be shortened if required. To do this, cut the protruding section. If necessary, the strap can be passed from the wrist to the hand.
- 6 Grasp the shoulder cap and place on the luxated clavicle (make sure that it is not positioned too far laterally) (fig. H). Guide the strap with the orange-coloured hook closure under the contralateral armpit and hook to the pin on the shoulder cap (fig. I). Close the relevant Velcro® closure; (fig. K); the rubberised end piece must be located precisely over the orange-coloured marking (this acts as a reference for the patient).
- 8 Tighten the contralateral straps with the Y-shaped Velcro® closures (fig. L, M, N). If the straps are too long, remove the Y-shaped Velcro® closure, cut the straps to length and replace the Y-shaped Velcro® closure again.
- 9 Tension the vertical strap on the back (fig. O), then the vertical strap on the front (fig. P).
- 10 Check to make sure all straps are tensioned.

REMOVING THE ORTHOSIS

- 1 Open the tensioning strap (this can be recognised by the large rubberised end piece; number [4]), but do not pull it out (fig. R).
- 2 Unhook the orange-coloured closure (strap number [3]) from the pin on the shoulder cap (fig. S) (if need be, first loosen the Velcro® closure with the rubberised end piece).
- 3 Open the wrist strap (fig. T).
- 4 Open the white clip buckle closure (fig. U) (waist strap) and carefully remove the orthosis; avoid sudden movements (fig. V).

PUTTING ON THE APPLIANCE FOR THE PATIENT

- 1 Make sure that the tensioning strap (this can be recognised by the large rubberised end piece; number [4]) is open.
- 2 Unhook only the white clip buckle closure (fig. A) (stomach strap, number [1]) and the orange-coloured hook closure (contralateral strap on the front, number [3]) (fig. B).
- 3 Position the arm pocket with forearm support so that it is at the height of the elbow (make sure it is not too low). The arm pocket must be on the outside with the Velcro® closure facing upwards.
- 4 Close the waist strap with the white clip buckle closure (fig. C).
- 5 Align the elbow pocket in such a way that it forms a line with the luxated collarbone; then guide the elbow into the arm pocket (fig. E).
- 6 Close the wrist strap (number [2]) around the wrist (fig. G).
- 7 Grasp the shoulder cap and place on the luxated collarbone (make sure that it is not positioned too far laterally) (fig. H). Guide the strap with the orange-coloured buckle under the contralateral armpit and hook into the pin on the shoulder cap (fig. I). Close the appropriate Velcro® closure (fig. K); tighten until the rubberised end piece is positioned over the orange-coloured marking.
- 9 Tighten the tensioning strap (number [4]) and secure with the Velcro® closure (fig. Q).

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