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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM**

Ankle support in AirX™ fabric

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/ EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.


WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Code	REF.A1181				
Size	S	M	L	XL	
Circ. of ankle from base of heel across the malleoli cm	26/29	29/32	32/35	35/38	
Colour	black				
fits right and left					

MAINTENANCE

 Do not bleach No chemical cleaning
 Do not iron Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Conservative treatment after distortions grade II
- Acute tendinopathies
- Postoperative rehabilitation
- Supports the resumption of sporting activities after distortions grade II and III
- Chronic joint pain caused by ankle osteoarthritis

CONTRAINDICATIONS

Currently no known

FEATURES AND MATERIALS

- Structure in breathable AirX™ fabric
- Upper closure system with stabilising straps and buckle
- Spiral reinforcing rods on both sides
- The tension straps are fixed crosswise across the arch and guided from underneath the sole of the foot to achieve maximum stability
- Opening at heel

PUTTING ON THE APPLIANCE

- 1 Open all Velcro® closures and temporarily close them slightly.
- 2 Pull the bandage over the foot so that the heel is exactly in the middle of the heel opening, and subsequently close the Velcro® closure strap (fig. A).
- 3 Positioning of the Velcro® closure on the inside:
 - grip the sewn Velcro® strap medially, guide over the instep and pull downwards to the opposite side (fig. B);
 - guide under the sole of the foot before the heel, tension and close vertically above the inner ankle (fig. C and D).
- 4 Positioning of the lateral Velcro® closure:
 - grip the sewn Velcro® strap laterally, guide over the instep and pull downwards to the opposite side (fig. E);
 - guide under the sole of the foot before the heel and close vertically above the outer ankle (fig. F and G).
- 5 Guide both Y-shaped ends of the fixing strap in the appropriate buckle, pull evenly with both hands and fix with the Velcro® trap (fig. H).

NOTE: The Velcro® straps can be tensioned to different extents, depending on which degree of stabilisation is required.

ONLY FOR THE FIRST FITTING

- 1 Shorten the Velcro® straps if necessary:
- 2 Remove the Y-shaped Velcro® straps at the end of the belt (fig. I). Cut the protruding part of the belt (fig. II).
- 3 Re-attach the Y-shaped Velcro® closure (fig. I).

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