
REF. 1336 W - REF. 1335 M

REF. 1636C W - REF. 1635H M - REF. 1635C M - REF. 1635H M

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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Padded orthosis for stabilising and supporting the dorsal-lumbar spine, with malleable reinforced stays

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly used.

SELECTION/SIZE

Item	REF. 1336 W - (REF. 1335 M)					
Size	XS	S	M	L	XL	XXL
Pelvic circumf. cm	75/85 (70/80)	85/95 (80/90)	95/105 (90/100)	105/115 (100/110)	115/125 (110/120)	120/135 (120/130)
Front height cm	21 (23)	21 (23)	21 (23)	21 (23)	21 (23)	21 (23)
Back height cm	55	55	55	55	55	55
Colour	beige					

Item	REF. 1636C W - REF. 1636H W (REF. 1635C M - REF. 1635H M)			
Size	S	M	L	XL
Pelvic circumf. cm	75/90 (70/85)	90/105 (85/100)	105/120 (100/115)	120/135 (115/130)
Front height cm	21 (23)	21 (23)	21 (23)	21 (23)
Back height cm	55	55	55	55
Colour	C: cream - H: blue			

Option: REF.9ST012 additional reinforced stay cm 38

MAINTENANCE
☒ Do not bleach ☒ No chemical cleaning
☒ Do not iron ☒ Do not tumble-dry
☒ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Dispose of this product and its parts responsibly

INDICATIONS
• Back pain due to arthritis or kyphosis
• After-effects of fractures of dorsal vertebrae and/ or upper lumbar vertebrae
• After-effects of stable compression fractures of the dorsal spine due to osteoporosis or metastasis

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS
• Traction system with elastic straps and loops, for targeted pressure and perfect stabilisation
• 4 linear back pockets allow the insertion of optional rigid reinforced stays as a replacement of the ones provided.
• Shoulder pads with axillary padding
• Velcro® front fastening
• Fabric made of cotton with a small amount of elastane
• Tissue in cotton containing a small quantity of elasthanne
REF. 1635 - 1636:
• New posterior insert of reticular fabric for maximum permeability

PUTTING ON THE APPLIANCE
1 Open up the orthosis and provisionally affix the Velcro® straps to the strips on the side (fig. A). Position the support in such a way that the inside label faces upwards.
2 Apply the support and use the middle Velcro® closure to fix in place securely (fig. B). Ensure that the closure is precisely centered.
3 Slowly tighten up the support by securing the lower Velcro® straps to the front strips (fig. C).
4 Then secure the upper straps to the front Velcro® strips (fig. D).
5 Place the upper tensioning straps over the shoulders, then feed under the arms and cross over at the back. Make sure they go through the loops (fig. E). Then secure to the front strips (fig. F).

To ensure optimal stabilisation, mould the rigid stays provided and slide them into the pockets provided in the back.

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