SaW 60°

REF. 84503

Schulterabduktionskissen

Shoulder abduction wedge

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Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares with sole responsibility that this product is a Class I medical product and has been manufactured in accordance with the basic requirements of the directives 93/42/EEC and 2007/47/EEC. These operating instructions were created in accordance with the specification of the above-mentioned directives. They aim to ensure the safe and correct use of the medical product.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SELECTION/SIZE

Size | universal |

Colour: black - fits right and left

Option: REF. 84505/P with pump

Structure of the cushion:
Cylinder-shaped, inflatable inner cushion made from robust plastic

Abduction of
45°/60°
[depending on
inflation level]

SCOPE OF DELIVERY

English

- A Adhesive Velcro® straps with closures and three-cornered eyelets
- B Inflatable inner cushion (optionally available with pump REF. 84505/P)
- C Cover with buckles and Velcro® closure
- D Wedge-shaped arm support and hand support
- E Abdominal strap
- F Lower arm securing straps

CHARACTERISTICS AND MATERIALS

- Structure of the cushion: Cylinder-shaped, inflatable inner cushion made from robust plastic
- Exterior made using gauze fabric with adhesive surfaces for the Velcro® closure
- Shoulder and pelvic strap made from nylon gauze with Velcro® closure
- Pad with sandwich structure with interior padding
- Wedge-shaped axilla support made from coated foam for improved stability and positioning of the cushion
- Abduction can be regulated by modifying the inflation of the cushion
- Lower arm and hand support

INDICATIONS

- Tuberculum-majus fracture with minor dislocation
- Treatment of capsulitis adhaesiva (frozen shoulder) after mobilisation
- Post-operative care of endoprosthesis and arthroprosthesis of the shoulder
- Post-operative care after surgery to the rotator cuff
- Nerve damages

CONTRAINDICATIONS

- Currently no known

ADAPTING FOR THE DOCTOR / TECHNICIAN

English

- 1 Insert the inner cushion into the cover which is usable with Velcro® (if the adhesive surface is already attached, remove it). Inflate the inner cushion to the required volume. The cushion must be positioned against the side of the patient to better ensure stability. Use the appropriate adhesive surface to close the cover (and conceal the valve).
- 2 Insert the arm of the healthy side through the ring formed by the shoulder and axilla strap so that the eyelet that is affixed to the adjustable strap faces forwards after the arm has been moved to the correct position (fig. A).
- 3 The eyelets sewn to the straps must, on the other hand, be on the back (fig. B).
- 4 Secure the cushion with the appropriate buckles to the straps (fig. C). Regulate the length of the straps so that the arm lies on the cushion with the desired angle of abduction of the shoulder.
- 5 Position the hand support in the most convenient holding position and secure it to the cover (fig. D). Provisionally attach both straps for the lower arm.
- 6 Grasp the hand support and secure the two lower arm straps (fig. E) in accordance with the doctor's instructions.
Note: Ensure that the head of the humerus is not pushed upwards.
- 7 Affix the wedge-shaped arm support under the medial (inner) side of the arm. (There should be no body contact, if possible.) Secure the wedge using the appropriate strap (fig. F).

PUTTING ON THE APPLIANCE FOR THE PATIENT

English

- 1 Insert the arm of the healthy shoulder through the ring formed by the shoulder and axilla strap so that the eyelet that is affixed to the adjustable strap faces forwards after the arm has been moved to the correct position (fig. A).
- 2 The eyelet sewn to the straps must, on the other hand, be on the back (fig. B).
- 3 Secure the cushion with the appropriate buckles to the straps (fig. C).
- 4 Grasp the hand support and secure the two lower arm straps (fig. D).
Note: Ensure that the head of the humerus is not pushed upwards.
- 5 Affix the wedge-shaped arm support under the medial (inner) side of the arm. (There should be no body contact, if possible.) Secure the wedge using the appropriate strap (fig. E).

MAINTENANCE

do not bleach

do not iron

Washing instructions: hand wash at max. 30° C with neutral soap.

Dry away from heat sources.

no chemical cleaning

do not tumble-dry

- Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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Headquarter: ORTHOSERVICE AG
Niederlassung Deutschland: Orthoservice Deutschland GmbH
Orthoservice AG (CH) · 6362 Stansstad (NW) · Switzerland