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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Back orthopedic brace

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. Please read the composition of the product on the inside label carefully. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

WARNINGS

The product has been designed for the specific indications stated below and must be prescribed by a doctor and fitted by an orthopaedic technician. The orthopaedic technician is responsible for fitting the product and providing information about how to use it safely in accordance with individual requirements. The product must be fitted with the utmost care to ensure its efficacy, tolerability and correct function. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/ orthopaedic technician be altered. The manufacturer is no longer liable if the product is used or adapted inappropriately. The orthosis is intended to be used only by a single patient; if it is used for more than one patient, the manufacturer accepts no responsibility or liability in accordance with medical device legislation. Direct contact with the skin may cause reddening or irritation in hypersensitive individuals. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician. The orthopaedic efficacy of the product may only be guaranteed when all of its components are used.

SELECTION/SIZE

Item	REF.91111
One size	Length 30 cm

MAINTENANCE

⚠ Do not bleach ⚠ No chemical cleaning
⚠ Do not iron ⚠ Do not tumble-dry
⚠ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Helpful for correcting a hunchback, postural kyphosis and scapula winging.
- A useful aid for back pain caused by poor posture.

CONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- The support and straps are made from PE foam bonded with velveteen and a technical fabric.
- A central loop at the back provides a better grip when fastening.
- The straps can be adjusted and shortened using the Y-Velcro® straps:
 - The overall length can be adjusted;
 - The positioned of the shoulder straps can be adjusted.
- Fastens with a snap-fit buckle.
- Easy to put on, just like a backpack.

PUTTING ON THE APPLIANCE

WEARING THE BRACE FOR THE FIRST TIME

- To put on the brace slide the patient's arms through the shoulder straps (fig A).
- Adjust the position of the shoulder straps using the Y-tab on the top end (fig B):
 - Position the shoulder straps so that the padding protects the armpit.
 - If they are too long, remove the Y-Velcro® strap, cut it to the correct length and reposition the Y-Velcro® strap.
- Adjust the snap-fit buckle at the front (fig C).
- Adjust the length of the lower straps:
 - Stretch the straps to fit the patient (fig D).
 - If they are too long, remove the Y-Velcro® strap, cut it to the correct length and reposition the Y-Velcro® strap (fig H + fig I).
- Place the brace back on the patient and adjust the fit if necessary.

WEARING THE BRACE ONCE FITTED

- Unfasten the snap-fit buckle.
- To put on the brace slide the patient's arms through the shoulder straps (fig E).
- Pull the straps at the sides (fig F) and fasten them at the front with the snap-fit buckle (fig G).

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