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PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM

Knee brace

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been prepared in its application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagneticfields (only REF. 33104, REF. 33116). Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Code	REF. 33160	REF. 33163	REF. 33104	REF. 33116		
Size <td>XS</td> <td>S</td> <td>M</td> <td>L</td> <td>XL</td> <td>XXL</td>	XS	S	M	L	XL	XXL
Thigh circ. taken 15 cm above mid-patella cm	33/39	39/45	45/52	52/58	58/66	66/75
Length cm GenuSKILL 04/16	34	34	34	34	34	34
Length cm GenuSKILL 60/63	27	27	27	27	27	27
Colour	grey					

GenusKILL 04/60/63: fits right and left - **GenusKILL 16:** specify right or left

MAINTENANCE

⚠ Do not bleach No chemical cleaning

🌀 Do not iron Do not tumble-dry

🧼 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Dispose of this product and its parts responsibly

INDICATIONS GENUKILL 04/60/63

- Knee pain due to bruises
- Early stage arthrosis
- Solo GenuSKILL 04/63: Chondropathy/ Chondromalacia patella
- Solo GenuSKILL 04: Mild patella instability
- Solo GenuSKILL 04/63: Prevention of patella luxation during sports
- Solo GenuSKILL 60: 1st grade sprains to collateral ligaments
- Solo GenuSKILL 60: Tendinitis

INDICATIONS GENUKILL 16

- Mid-lateral patella instability
- Patellofemoral arthrosis
- Gonarthrosis

KONTRAINDICATIONS

Currently no known

CHARACTERISTICS AND MATERIALS

- Structure made of BreathSkill:
 - soft and breathable
 - 3D weave
 - high resistant mesh
- CompresSkill back support:
 - stretchy multi-fiber and compression material
 - low profile anti-irritation (knee flexion)
 - fully adaptable
 - maximum adherence
 - parfaite adhérence au membre
- No open kneecap GenuSKILL 60 only
- Reduced height GenuSKILL 60 / 63 only
- Open kneecap GenuSKILL 63 only
- FlexSkill front support (bi-layer technical fabric) with special low profile anti-irritation welding seam
- Patella stabilizer included GenuSKILL 04 only
- Medial and lateral stays to keep brace correctly positioned. GenuSKILL 04 /16 only
- Patella alignment system GenuSKILL 16 only
 - lateral J-shaped foam stabilizer
 - patella support and pressure system with winches (adjustable)

PUTTING ON THE APPLIANCE

GenusKILL 60

1 Slide the tubular support onto the knee with the orange label towards the top (fig. A).

2 Position the knee support so that the CompresSkill fabric is behind the knee.

GenusKILL 04 - GenuSKILL 63

1 Position the tubular support onto the knee with the orange label towards the top (fig. B).

2 Position the knee support so that the kneecap is perfectly centred in the opening provided (fig. C).

GenusKILL 16

1 Undo the patella straps and fasten them back on themselves for the time being.

2 Slide the tubular support onto the knee with the orange label towards the top (fig. D) taking care to align the opening with the kneecap. Lateral instability: the stabiliser must be lateral to the patella. Medial instability: the stabiliser must be medial to the patella (fig. E).

3 Pass the upper strap through the buckle provided, tighten, and fasten with the Velcro® (fig. F).

4 Repeat this procedure with the lower strap (fig. G).
N.B.: if the straps are too long, they can be cut to size and the corresponding Y-shaped Velcro® can be reattached (fig. H, I).

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Velcro® è un marchio di fabbrica registrato di Velcro Industries B.V.