# regain

REF. 95050 (corta)  
REF. 95000 (lunga)

Postoperative Kniegelenkorthese  
Post-operative knee brace  
Genouillère postopératoire  
Послеоперационный наколенник  
Ginocchiera post-operatoria a movimento graduato

## DEUTSCH

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## ENGLISH

Please read the instructions carefully and thoroughly

### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

### INDICATIONS

- Postoperative care when protection of the flexion and extension of the knee is required (ligament reconstruction, transposition of the anterior tibial apophysis, fracture of the distal third of the femur, fractures of the proximal third of the tibia, fractures of the patella, tendon injuries of the extensor system)
- Post-traumatic care of knee distorsions
- Valgisation and varisation osteotomies

### CONTRAINDICATIONS

Currently no known

### MAINTENANCE

- Do not bleach  No chemical cleaning
- Do not iron  Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap, dry away from heat sources.

#### Washing instructions

- Frame: Clean using a sponge with lukewarm water and a neutral soap Then dry with a cloth
- Padding: Remove the padding and clean in lukewarm water with a neutral soap. Dry away from heat sources.
- Check screws for tight fit regularly
- Do not dispose of the product or any of its components into the environment.
- Contact an orthopaedic technician for the replacement of used parts

### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

### PRECAUTIONS

It is advised that the pressure exerted by the device does not act on parts of the body suffering from wounds, swelling or lumps. It is advisable to wear a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopedic technician. The device should not be worn in the vicinity of open flames or strong electromagnetic fields.

### WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopaedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopaedic technician. Under no circumstances should the adjustment made by the physician/orthopaedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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## ENGLISH

### CHARACTERISTICS AND MATERIALS

- Structure with easily malleable, robust and anatomically adaptable hinged bars made from light aluminium alloy
- MTP pads prevent slipping of the orthosis. Only covered on the outside with non-toxic, latex-free and hypoallergenic material
- Protective condyloid pads, removable and washable
- Monocentric steel joint
- New system for setting flexion and extension PULL&TURN; it permits easy setting with minimum effort.
- Extension: 0° to 70° in 10°-increments - flexion: 0° to 120° in 10°-increments
- Closure system with Velcro® straps ensuring a safe fit of the orthosis on the leg
- Possibility of shortening the orthosis by means of guide marks on the bars
- Pads can easily be cut to size for better adaptation to the circumference of the leg

### ADAPTING

#### ADAPTING FOR THE DOCTOR / TECHNICIAN

- 1 Remove the pads from the splints with the belts.
- 2 Attach the pads in such a way that they wrap tightly around the leg, and close at the front using the appropriate Velcro® closures. First close the ones at the thigh, then at the calf (fig. A). If necessary, cut the pads to fit the patient's leg (only Regain II: leave out the knee /only Regain II open: cut along the pre-pressed cutting lines so that the padding does not open).
- 3 Mould the splints (by bending) and adjust to the leg profile of the patient (fig. B).
- 4 Carry out the flexion and extension settings as described below (points 9-15): Carry out the same procedure for both joints.
- 5 Fix the splints at the pads (fig. C): The splints have to run along the centre line of the leg, the joints have to be at the height of the patella so that they rest on the condyles. Use the lateral ankle and the trochanter major as reference points when aligning the splint with the centre line of the leg.
- 6 Only Regain II open: Ensure that the straps are positioned exactly on the pads.
- 7 Ensure that both joints are at the same height.
- 8 Guide the Velcro® straps through the appropriate rings and close them (fig. D): First fasten the straps of the thigh. Start with the straps closest to the knee, then proceed to fastening the straps of the calf. Tighten the Velcro® straps so that the orthosis cannot slip vertically.

**NOTE:** Bending has to be done manually. A bending iron can be used, or the orthosis can be put over a rounded object (chair, edge of the bed...).

#### SETTING THE JOINTS (for the physician/orthopaedic technician)

- 9 Pull out the locking mechanism (A) by pressing the push buttons on the side and turning the cap (B) simultaneously. Release the locking mechanism (A).
- 10 Turn the cap (B) until the opening is on top of the flexion stop (C).
- 11 Pull the flexion lock (C) outwards.
- 12 Turn the cap (B) and move the flexion lock until it reaches the desired degree of flexion (the degrees of flexion are indicated on the mechanism).
- 13 Push the flexion lock (C) inwards until it locks in place in the selected position.
- 14 Repeat steps 10 to 13 also for setting the extension by using the lock (D).
- 15 Turn the cap (B) until the safety lock engages (A).

#### RECOMMENDATIONS FOR THE PATIENT

- To facilitate washing the immobilized limb without removing the brace, proceed as follows:
- 1 Open the upper part (straps and padding) and wash the thigh
  - 2 Close the upper part
  - 3 Open the lower part (straps and padding) and wash the leg
  - 4 Close the lower part
  - 5 Tighten the strap to prevent vertical migration of the brace.
- To wash the limb we recommend using a damp sponge.

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