podifer 50

DÄMPFUNGSEINLAGE MIT QUERVERSTÄRKUNG
CUSHIONING INSOLE WITH TRANSVERSE REINFORCEMENT
SEMELLE AMORTISSEUSE À RENFORT TRANSVERSAL
АМОРТИЗИРУЮЩИЙ СУПИНАТОР

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Modern ultra-thin insoles for people that spend long hours walking or standing and want more relaxed feet and enhanced comfort when walking. The special shock-absorbing material in the heel area cushions the energy produced as the foot touches the floor. An optimum combination of foam layers ensures the same effect during the rest of the stance phase. A PodoDry layer wicks sweat away from the foot, while the special PA and elastane covering extends the life of your insoles.

INDICATIONS: Inflammation of the heel * Plantar fascitis * Flatfoot or high arched foot * Metatarsalgia.

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INDICATIONS: Inflammation of the heel * Plantar fascitis * Flatfoot or high arched foot * Metatarsalgia.

Cosmopolitan modern ultra-thin insole for those who walk a lot or stand a lot and want more relaxed feet and enhanced comfort when walking. The special shock-absorbing material in the heel area cushions the energy produced as the foot touches the floor. An optimum combination of foam layers ensures the same effect during the rest of the stance phase. A PodoDry layer wicks sweat away from the foot, while the special PA and elastane covering extends the life of your insoles.

INDICATIONS: Inflammation of the heel * Plantar fascitis * Flatfoot or high arched foot * Metatarsalgia.

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ENG. PUTTING ON THE APPLIANCE

FIRST APPLICATION

- Remove the original insole from the patient's shoe.
- Place the insole in the shoe, check that it is correctly positioned, and ask the patient to try it on.
- Alterations may be made to the insole to suit individual requirements (it is advisable for the patient to wear the insole for a few hours before making any alterations). The front can be shaped with scissors (use the insole that was previously in the patient's shoe for reference). The arched detailing supporting the longitudinal arch can be filed down a little at a time, as can the outer edge of the rear part of the insole.
- We recommend removing the insole from the shoe after wear so that it can air.

SUBSEQUENT WEAR

- Place the insole in the shoe and check that it is correctly positioned.

- We recommend removing the insole from the shoe after wear so that it can dry.

CONTRAINDICATIONS

Currently no known

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

MANTENANCE

Washing instructions:

Wash by hand up to 30°C with neutral soap; let it dry away from heat sources. Do not dispose of the product or any of its components into the environment.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician.

Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

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