REF. 20414

Orthese zur Immobilisierung der Schulter
Orthosis for shoulder immobilization
Orthèse d'immobilisation des épaules
Ортез иммобилизации плеча
Ortesi per immobilizzazione spalla

Headquarter: ORTHOSERVICE AG
Via Milano 7 - CH-6830 Chiasso (TI) - Switzerland
Tel. 0041 (0) 91 822 00 88 - Fax 0041 (0) 91 822 00 89
info@orthoservice.com - www.orthoservice.com

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ENGLISH

Please read the instructions carefully and thoroughly

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG states under its sole responsibility that this product is a Class I medical device that has been manufactured in accordance with the requirements of Directives 93/42/EEC and 2007/47/EC. These instructions have been written in conformity with these directives, and their aim is to ensure the safe and proper use of this medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro® Industries B.V.

SELECTION/SIZE

Item: REF.20414	S	M	L	XL
Size	S	M	L	XL
Forearm length	28/34	34/42	42/48	48/58

Colour: black / Fits right and left

*measurement taken from the elbow to the end of the closed fist

MAINTENANCE

- Do not bleach
- No chemical cleaning
- Do not iron
- Do not tumble-dry
- Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

PRECAUTIONS

It is advised that the pressure exerted by the product does not act on parts of the body suffering from wounds, swelling or lumps. The product should not be tightened too firmly so as to avoid the development of pressure points or compression of the nerves or blood vessels underneath. It is advisable to wear the product over a garment, avoiding direct contact with the skin. If in doubt about how to apply, contact an orthopaedic technician. The device should not be worn in the vicinity of open flames.

WARNINGS

The product may only be used in the event of medical indication or when prescribed by a physician and must be applied by an orthopedic technician, the competent authority figure for the application as well as for information on the safe use, according to individual needs. To ensure the effectiveness, tolerability, and correct operation, it is necessary that the device is put on with the utmost care. Any changes to the structure or adjustment of the device must be prescribed by a physician and performed by an orthopedic technician. Under no circumstances should the adjustment made by the physician/orthopedic technician be altered. Use of the product is recommended for a single patient. With hypersensitive patients, direct contact with skin may cause redness or irritation. In the case of onset of pain, swelling, lumps or any unusual reaction, immediately contact your physician.

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ENGLISH

CHARACTERISTICS AND MATERIALS

- Poly cotton fabric sling with an elbow support.
- Immobilising chest strap with Velcro® fastening and padding in the central part.
- Shoulder strap with Velcro® fastening and padding in the central part.

INDICATIONS

- Simple humeral fractures
- Glenohumeral luxation
- Distorsions of the shoulder
- Postoperative care of humeral fractures and interventions due to dislocations of the shoulder
- Postoperative care for tendon sutures of the long head of the biceps muscle

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SUBSEQUENT WEAR

- Undo the Velcro® fastenings on the chest strap and shoulder strap.
- Place the bent arm in the fabric sling (fig. 10).
- Pass the shoulder strap (black) over the opposite shoulder (fig. 11a), place the tip through the corresponding black ring, and fasten with Velcro® (fig. 11b).
- Pass the (white) chest strap around the patient's chest, place the tip through the corresponding white ring (fig. 12a), and fasten with Velcro® (fig. 12b).

WARRANTY

This warranty is governed by the statutory provisions of the country in which the product was purchased. Where relevant that country's regulations concerning warranties and the relationship between the seller and buyer shall apply. If you believe that the terms of the warranty may be applicable, first contact the retailer where the product was purchased. Orthoservice will repair or replace all or part of the device and its accessories if defects in materials and workmanship arise during the period of six months from the date of purchase.

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