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 **PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM**

Knee orthosis for post-surgery with graduated movement

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is also advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels.

If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields. Do not apply in direct contact with open wounds.

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.G1042 • REF.G1043		REF.G1041
	Regular	Long	Rehab
Size	< 81	> 81	
Limb length groin - ankle cm			
Length cm	56	65	40
Colour	black		

Fits right and left

MAINTENANCE

-  Do not bleach
-  Do not iron
-  Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.
-  No chemical cleaning
-  Do not tumble-dry

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Post-surgery for the protection of knee flexion and extension (ligament reconstruction, transposition of the anterior tibial apophysis, fracture of the distal third of femur, fractures of proximal third of tibia, fractures of patella, tendon injuries of extensor system)
- Post-traumatic stage of knee sprains
- Valgization and varization osteotomies

CHARACTERISTICS AND MATERIALS

- Valves with strap-loops and structural hook-and-loopsystem
- Non-metallic adjusting joint, light and intuitively adjustable
- Flexion: from 0° to +120° with intervals of 10°
- Extension: from -10° to +60° with intervals of 10°
- Possibility to lock joint from 0° to +60° with intervals of 10°
- Removable and washable condyle protection pads
- Structure with hinges made of light aluminum alloy
- Possibility of shortening brace through marks on the hinges
- MTP Pads to limit brace migration and protect leg from direct contact with hinges

ADAPTING

- 1 Open all straps and temporarily close them upon themselves.
 - 2 Separate the padding from the rods with the straps.
 - 3 Apply the padding such that it fits closely to the leg, closing it at the back with the appropriate Velcro® fasteners: first the thigh padding, then the calf padding (fig. A). In order to position the padding at the optimal height, place the orthosis next to the leg. If necessary the padding may be shortened to adapt it to the patient's legs (only for G1043: shorten the padding along the thinner parts in such a way that the padding covers properly).
 - 4 Model the rods (by bending them) to fit them to the shape of the patients' legs (fig. B).
 - 5 Adjust the flexion extension as described (points 1-5, hinge adjustment): perform the same operation on both joints.
 - 6 Fix the rods to the pads (fig. C): the rods should follow the median line of the leg and the joints should be aligned with the patella so as to rest on the condyles. To align the rods with the median line of the leg, use the lateral malleolus and the large trochanter as reference points.
 - 7 Ensure that the Velcro® straps are perfectly aligned with the corresponding paddings.
 - 8 Ensure that the two joints are at the same height.
 - 9 Verify that the orthosis is correctly fitted in the appropriate rings (fig. D-E): fasten first the thigh straps beginning with the one closest to the knee; then proceed by fastening the calf straps. Tighten the straps so as to avoid vertical migration of the guide.
- NOTE:** The bending procedure should be performed by hand, with the help of the grips or laying the orthosis on a rounded body (chair, edge of a bed); do not use your own knee as a lever.

ADAPTING FOR THE DOCTOR / TECHNICIAN

HINGE ADJUSTMENT - REF. G1042/G1043

- 1 Press the orange button and at the same time rotate the transparent cover (fig. 1).
- 2 To adjust flexion (fig. 2):
 - Rotate the opening of the cover to coincide with the flexion stop and remove the pin
 - Rotate the cover and the pivot to the hole with the degree of flexion required (shown on the mechanism) and insert the pin
- 3 Rotate the cover again (fig. 3).
- 4 To adjust extension (fig. 4):
 - Rotate the opening of the cover to coincide with the extension stop and remove the pin
 - Rotate the cover and the pivot to the hole with the degree of extension required (shown on the mechanism) and insert the pin.
- 5 Rotate the cover until the orange safety button clicks (fig. 5).

ADAPTING FOR THE DOCTOR / TECHNICIAN

SHORTENING THE RODS (OPTIONAL) - REF. G1042/G1043

- If necessary shorten the rods (the guide markings are shown in fig. 6):
- 1 Free the rods from the padding (thigh or calf) (fig. 6).
 - 2 Remove the straps from the terminal and intermediate valves (fig. 7).
 - 3 Remove from the two rods the valves pressing towards the centre, placing the screwdriver in the appropriate notch and pressing outwards (fig. 8).
 - 4 Place the rod on a worktable so that the edge is aligned with the grooves on the rod (which coincides with the desired length). Press downwards till it breaks (fig. 9).
 - 5 Reposition the terminal valve on the shortened rod, which need not be smoothed down. Slide in the straps.

RECOMMENDATIONS FOR THE PATIENT

- To facilitate washing the immobilized limb without removing the brace, proceed as follows:
- 1 Open the upper part (straps and padding) and wash the thigh
 - 2 Close the upper part
 - 3 Open the lower part (straps and padding) and wash the leg
 - 4 Close the lower part
 - 5 Tighten the strap to prevent vertical migration of the brace. To wash the limb we recommend using a damp sponge.

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