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 PLEASE READ THESE INSTRUCTIONS CAREFULLY AND SAVE THEM

Foam cervical collar

DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames. Do not apply in direct contact with open wounds

WARNINGS

It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter the adjustment made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall lapse in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices.In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe adverse event, report the incident to the manufacturer and the competent authority of your country. The effectiveness of the orthopedic product is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF.8300 - REF.8301
Size	universal
Height cm	7,5 10
Colour	white

Item	REF.8302
Size	S M L XL
Neck circ. cm	30/35 36/41 42/47 48/53
Height cm	6-8-10 6-8-10 6-8-10 6-8-10
Colour	white

Item	REF.8304
Size	S S S S M M M M L L L L
Neck circ. cm	28/34 29/35 30/36 31/37 30/36 31/37 32/38 33/39 36/44 37/46 38/48 39/49
Length in cm	47 49 50 52 50 52 54 56 58 60 62 64
Height cm	6 8 10 12 6 8 10 12 6 8 10 12
Colour	blue

MAINTENANCE

⚠ Do not bleach ☒ No chemical cleaning

⚠ Do not iron ☒ Do not tumble-dry

⚠ Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

INDICATIONS

- Slight distortions of cervical spine (whiplash)
- Torticollis
- Cervical vertigo
- Arthritis of cervical spine

CONTRAINDICATIONS

- Currently no known

CHARACTERISTICS AND MATERIALS

- REF. 8304 only:** Semi-rigid foam structure (with rounded edges)
- Covered with tubular jersey
- Back Velcro® fastening system
- 100% washable (warm water (30°C) and neutral soap)
- Radiolucent
- REF. 8304 only:** Perfect anatomical shape
- REF. 8300 - 8301 only:** External band reinforcement made of non toxic leather-like material

PUTTING ON THE APPLIANCE

- Hold the cervical collar in such a way that the concave side points upwards and the belt is on the outside.
- Place the concave front part of the cervical collar under the chin.
- Close the cervical collar at the back using the imitation leather band.

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