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
PLEASE READ THESE INSTRUCTIONS CAREFULLY AND KEEP THEM


Brace for hip and coxalgia

DECLARATION OF CONFORMITY
As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

TRADEMARKS OF MATERIALS
Velcro® is a registered trademark of Velcro Industries B.V.

SAFETY PRECAUTIONS
We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts. It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment that avoids direct contact with the skin. In the case of any change in symptoms, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully. We recommend to not wear the device in the proximity of free flames or strong electromagnetic fields, or to apply in direct contact with open wounds.

WARNINGS
It is advisable that the device, designed for the specific indications as set out below, be prescribed by a doctor or a physiotherapist and applied by an orthopedic technician, according to the specific needs of the patient. To ensure effectiveness, tolerability, and proper functioning, the application must be carried out with the utmost care. Never alter adjustments made by the doctor/physiotherapist/ or orthopedic technician. The responsibilities of the manufacturer shall remain in case of inappropriate use or adaptation. The orthosis is made and to be used by only one patient. Should the orthosis be inappropriately used, the manufacturer declines all responsibility, as provided by the regulation for medical devices. In hypersensitive individuals, direct contact with the skin may cause redness and irritation. In case of pain, swelling, welts, or any other adverse reaction, contact your doctor immediately, and in case of a particularly severe allergic reaction to any of the products, contact the manufacturer and the competent authority of your country. The effectiveness of orthopedic products is only guaranteed when all the components are properly in use.

SELECTION/SIZE

Item	REF. 91006 (Hipocross)
Size	S M L XL
Pelvis circ. cm	72/93 93/116 116/140 140/165
Circ. waist cm	65/84 84/109 109/140 140/170
Proximal thigh circ. cm	40/54 50/65 60/76 70/87

Colour: black / specify right or left

SELECTION/SIZE

Item	REF. 91007 (Adjustable rod for Hipocross)
Size	universal

Colour: black / specify right or left

MAINTENANCE
 Do not bleach No chemical cleaning
 Do not iron Do not tumble-dry
 Washing instructions: hand wash at max. 30° C with neutral soap. Dry away from heat sources. Dispose of this product and its parts responsibly

INDICATIONS
- Post-arthrosopic care.
- Post-surgical care (compatible with the degree of stabilization required).
- Conservative treatment of hip dislocations.
- Arthritis of the hip joint.
- Coxalgia or inoperable joint misalignment.

CONTRAINDICATIONS
Currently no known

CHARACTERISTICS AND MATERIALS
- Made of breathable AirX™ fabric.
- Lateral and anterolateral spiral support stays.
- Front Velcro® fasteners.
- Colour-coded Velcro® areas for fastening the straps.
- Elastic calibrated rods, with lateral crossing, for:
- Prevention of excessive adduction
- Partial limitation of flexion-extension
- Gel insert to prevent oedema
- A rigid bar that locks the joint in flexion and neutral abduction during the initial days following a dislocation
- REF.91007 – Adjustable rod for Hipocross:
- restrict extension (0° to 120° in 10°-increments)
- restrict flexion (0° to 120° in 10°-increments)
- regulate adduction.

PUTTING ON THE APPLIANCE FIRST APPLICATION
1 Position the brace on the hip with the vertical Velcro® facing outwards along the median of the body (fig. A)
2 Fasten the Velcro® pelvis strap, checking that the fabric is stretched out properly (fig. B)
3 Fasten the Velcro® thigh strap; hold the loop at the end of the strap if this is easier and check that the fabric is stretched out properly (fig. C)
4 Attach the elastic straps on the lower parts of the pelvis (blue – fig. D) and thigh (orange – fig. E)
5 Fasten the Velcro® fasteners on the sides (on the trochanters) and fasten the strap with the orange tip at the front on the pelvis (opposite iliac crest) (fig. F); fasten the strap with the blue tip on the thigh (fig. G). The fit can be tightened or loosened using the hooks at the front and back.
6 Use the bar provided to lock the flexion and abduction of the joint in the neutral position; attach the bar to the side strip (fig. H) and fasten the straps provided around the thigh and pelvis (fig. H)


REF.91007
7 To use the adjustable rod REF.91007, attach it with Velcro® to the side strip (fig. I)
8 Fasten the straps around the thigh and pelvis (fig. L)
9 The elastic straps can also be fastened over the stay as described above for a tighter fit (fig. N)

ADJUSTMENT OF THE JOINTS
10 Extension adjustment (fig. O)
11 Flexion adjustment (fig. P)
12 This washer (fig. Q) can be turned so that the respective wedge (fig. R) locks in the adjusted flexion-extension setting (fig. S) for the screw provided (to hold the wedge in place)

SETTING THE ADDUCTION/ABDUCTION
14 Insert the supplied key into the appropriate seat and rotate until the desired abduction / adduction degree is achieved (fig. T)

SUCCESSIVE APPLICATIONS
1 Position the brace on the hip with the vertical Velcro® facing outwards along the median of the body (fig. 1)
2 Fasten the Velcro® pelvis strap, checking that the fabric is stretched out properly (fig. 2)
3 Fasten the Velcro® thigh strap; hold the loop at the end of the strap if this is easier and check that the fabric is stretched out properly (fig. 3)
4 Hold the elastic straps, cross them at the sides and fasten: the strap with the orange tip at the front on the pelvis (opposite iliac crest) (fig. 4); fasten the strap with the blue tip on the thigh (fig. 5).

COOLING
1 Take the gel cushion out of its velour cover
2 Place the cushion into the freezer at least 2-3 hours prior to its application (note: do NOT put the velour cover in the freezer)
3 While applying, insert the gel cushion in the velour cover.
4 Insert the gel cushion and the velour cover in the special pocket inside the brace (fig. 6).
5 Apply the brace

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