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Retain the configuration as adjusted by the doctor/technician. In order to remove the orthotic, open the arm fixation and unhook the push-in buckle fasteners. To put the orthotic on again, proceed as follows:

- Position the pillow on the hip on the same side as the shoulder requiring immobilisation and rest the arm in the sling.
- Fasten the waist strap with the side-release buckle (fig.B).
- Close the sling with the Velcro fasteners and slide the thumb into its support.
- Pass the vertical strap over the opposite shoulder and fasten the siderelease buckle (fig.E).

#### PUTTING ON THE APPLIANCE FOR THE DOCTOR/TECHNICIAN

1 Undo the two side-release buckles. (REF. 410 ONLY) If necessary, detach the sling from the pillow and reposition it on the actual pillow so it matches up with the marking for the desired degree of abduction: orange=10°; white=15° (fig. A1). (EXCEPT FOR REF. 410) If necessary, detach the sling from the pillow and reposition it with Velcro® on the actual pillow on the same side as the shoulder requiring immobilisation. (fig. A2)

2 Position the pillow on the hip on the same side as the shoulder requiring immobilisation.

3 Fasten the waist strap with the side-release buckle (fig. B) and tighten with the corresponding Velcro® (fig. C).

4 Carefully slide the forearm into the sling (fig. D1), fasten it with the Velcro® strap provided and (REF. 410 ONLY) adjust the length of the back strap (fig. D2).

### Shoulder brace with immobilizer in fixed abduction

#### DECLARATION OF CONFORMITY

As the manufacturer, ORTHOSERVICE AG declares, under its sole responsibility, that this is a class I medical device that has been manufactured and assessed, according to the EU Regulation 2017/745 (MDR). These instructions have been drawn up in application of the aforementioned Regulation. They are intended to ensure an adequate and safe use of the medical device.

#### TRADEMARKS OF MATERIALS

Velcro® is a registered trademark of Velcro Industries B.V.

#### SAFETY PRECAUTIONS

We recommend that the pressure exerted by the device does not act on parts of the body with wounds, swelling, or welts.

It is advisable to not over tighten the device so as to avoid generating areas of excessive local pressure or the compression of underlying nerves and/or blood vessels. Furthermore, we recommend wearing a garment to avoid direct contact with the skin. If in doubt on how to apply the device, contact a doctor, physiotherapist, or an orthopedic technician. Read the product composition on the internal label carefully.

We recommend to not wear the device in the proximity of flames. Do not apply in direct contact with open wounds.

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Colour: black - fits right and left

\*Measurement taken from elbow to end of closed fist

#### Code

Code	REF. E110 10° - REF. E210 15°	REF. E111 10° - REF. E211 15°
Size	universal	universal

Colour: black - specify right or left

#### MAINTENANCE

- ✗ Do not bleach
- ✗ No chemical cleaning
- ✗ Do not iron
- ✗ Do not tumble-dry

✓ Waschanweisung: Wash the cushion cover (after the foam inner padding has been removed), belts and the neck and underarm padding with hand with warm water (at max. 30°C) and neutral soap, rinse thoroughly. Dry away from heat sources.

Do not dispose of the product or any of its components into the environment.

#### INDICATIONS

- Post surgery (in abduction) to repair rotator cuff and soft tissue, glenoid labrum suture, capsular stabilization

#### CONTRAINDICATIONS

Currently no known

#### CHARACTERISTICS AND MATERIALS

- Pad: high density foam with silver fiber in contact with the body
- Pad with front edge to improve stability;
- Only REF.410 Tilted toward arm support for a better comfort
- Protection neck sleeve
- External gauze nylon cover and sleeves with Velcro® fastening system
- Breathable arm pocket material
- Fixing straps with adjustable buckles with Velcro® fastening system
- Supporting strap for thumb

#### PUTTING ON THE APPLIANCE

##### ADAPTING FOR THE DOCTOR/TECHNICIAN

1 Undo the two side-release buckles. (REF. 410 ONLY) If necessary, detach the sling from the pillow and reposition it on the actual pillow so it matches up with the marking for the desired degree of abduction: orange=10°; white=15° (fig. A1). (EXCEPT FOR REF. 410) If necessary, detach the sling from the pillow and reposition it with Velcro® on the actual pillow on the same side as the shoulder requiring immobilisation. (fig. A2)

2 Position the pillow on the hip on the same side as the shoulder requiring immobilisation.

3 Fasten the waist strap with the side-release buckle (fig. B) and tighten with the corresponding Velcro® (fig. C).

4 Carefully slide the forearm into the sling (fig. D1), fasten it with the Velcro® strap provided and (REF. 410 ONLY) adjust the length of the back strap (fig. D2).

5 Pass the vertical strap over the opposite shoulder and fasten the siderelease buckle (fig. E).

6 Tighten the shoulder strap with the corresponding Velcro®. If it is too long, cut to size and reposition the Y-shaped Velcro® strap.

7 Position the padded sleeve to protect the neck.

8 Position the thumb support (and exercise ball) and attach it to the pillow with the Velcro® (fig. F).

#### PUTTING ON THE APPLIANCE FOR THE PATIENT

Retain the configuration as adjusted by the doctor/technician. In order to remove the orthotic, open the arm fixation and unhook the push-in buckle fasteners. To put the orthotic on again, proceed as follows:

- Position the pillow on the hip on the same side as the shoulder requiring immobilisation and rest the arm in the sling.

2 Fasten the waist strap with the side-release buckle (fig.B).

3 Close the sling with the Velcro fasteners and slide the thumb into its support.

4 Pass the vertical strap over the opposite shoulder and fasten the siderelease buckle (fig.E).

5 Tighten the shoulder strap with the corresponding Velcro®. If it is too long, cut to size and reposition the Y-shaped Velcro® strap.

6 Position the padded sleeve to protect the neck.

7 Position the thumb support (and exercise ball) and attach it to the pillow with the Velcro® (fig. F).

8 Position the shoulder strap with the corresponding Velcro® and attach it to the pillow with the Velcro® (fig. E).

9 Position the vertical strap with the corresponding Velcro® and attach it to the shoulder strap with the Velcro® (fig. D).